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Chargennummer

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Hersteller

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Instructions for Use partial denture alloy

ADORBOND® CC PLUS Powder 10 µm – 30 µm

ADORBOND® CC PLUS Powder 10-30 µm is a dental partial denture alloy based on cobalt. ADORBOND® CC PLUS Powder 10-30 µm is free from nickel, cadmium, beryllium and lead and fulfils the standards of EN ISO 22674 type 5 for appliances in which parts require the combination of high stiffness and strength, e.g. thin removable partial dentures, parts with thin cross-sections and clasps. ADORBOND® CC PLUS Powder 10-30 µm can be used as a partial denture alloy as well as a metal-ceramic alloy.

Composition w_i

Co	%	63,6
Cr	%	24,8
W	%	5,5
Mo	%	5,0
Si	%	1,1

Properties

Particle size d ₁	µm	10-30
Density ρ	g · cm ⁻³	8,2
Vickers hardness	HV 10	520
Melting range T _S - T _L	°C	1360-1420
0,2-% Yield strength R _{p0,2}	MPa	1020
Modulus of elasticity E	GPa	275
Tensile elongation at break A ₅	%	5

Recommendations for Use

Digital Design Modellation

The design is carried out with suitable CAD software in consideration of dental rules. Wall thickness after finishing of the object: at least 0,5 mm. Base thickness of standard bases: 0,75 – 0,9 mm and 1,0 – 1,2 mm for smaller bands or skeletonised frames. Thickness of lower jaw brackets: 1,8 – 2,0 mm x 4,0 x 4,2 mm. Thickness of transition area from minor connectors to clasps: 1,5 x 2,0 mm. Thickness of clasp tip: 1,2 – 1,5 mm. Round off the areas of transition from minor connectors to clasps and the tips of clasps.

SLM-Process

Avoid dust production while opening the powder container, transportation and filling into the SLM system! Use a suitable SLM system for the particle distribution of ADORBOND® CC PLUS Powder and carry out the SLM-Process according to the manufacturer's instructions for the SLM system. Non melted powder can be reused but should be sieved with a powder sieve (80 µm) first. Melted powder cannot be reused by e.g. casting.

Shot Peening

Blast the production platform with spherical ceramic material with Ø 0,125-0,250 mm and 2,5-3,5 bar until a uniform saturation of the surface is visibly achieved.

Stress Relief Heat Treatment

After the laser sintering process has been carried out, the building panels must be heat-treated, which reduces internal stresses. A suitable furnace with inert gas (argon) or vacuum function is recommended for this purpose. The heat treatment needs to be carried out in a special furnace under argon atmosphere. Place the production platform in an inert gas box with an argon flow of 1-2 L/min and put the box into the cold furnace. Then:

1. Within 60 min heat up to 450 °C and hold the temperature for 45 min.
2. Within 45 min heat up to 750 °C and hold the temperature for 40-60 min
3. Let the platform cool down. Open furnace door at max. 100 °C and remove the inert gas box and stop the argon flow.

Note: It is possible to carry out the test without inert gas, in which case increased attention must be paid to the fact that increased oxide formation may occur. The cooling phase should be carried out in the oven until room temperature is reached.

Treatment of the restoration

After the heat treatment and cooling down:

1. Remove restorations from the plate using a band saw, rotating instruments or pliers. Use pliers to remove the rest of the supports. Carbide cutters are recommended for finishing of the object.
2. Sand blast the surface by use of a pencil-blaster with aluminium oxide 100 µm or 250 µm.
3. Electrolytic polishing with commercial electrolytes in dental polishing units. Cover clasps and fitting parts with covering varnish while polishing.

Finishing

After finishing and fitting smooth the frame with a rubber polisher.

Soldering and Welding

Soldering of the frame can be carried out with ADORBOND® CC-Lot and high temperature flux. The width of the solder gap should be 0,05-0,2 mm. After brazing, flux residues and metal oxides must be acidified and the surfaces cleaned by steam blasting. For welding with a laser use ADORBOND® CC-Laserdraht.

Safety Note

Metal dust is harmful to health and can cause irritation when inhaled and in contact with skin. When handling the powder, working out and blasting, use suction and wear protective goggles, protective gloves, respiratory mask (according to FFP3 - DIN EN 149) and protective clothing. After working with the powder or the restorations, thoroughly clean the skin contact areas, hands are particularly important.

Warranty

These application recommendations are based on own experiments and experiences and can therefore only be regarded as guidelines. The dentist or dental technician is responsible for the correct processing of this alloy.

Batch number

Refer to instructions for use

Manufacturer

Not for reuse

Instructions for Use metal-ceramic alloy

ADORBOND® CC PLUS Powder 10 µm – 30 µm

ADORBOND® CC PLUS Powder 10-30 µm is a dental metal-ceramic alloy based on cobalt. ADORBOND® CC PLUS Powder 10-30 µm is free from nickel, cadmium, beryllium and lead and fulfils the standards of EN ISO 22674 type 5 for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered crowns, wide-span bridges or bridges with small cross sections, bars, attachments and implant retained superstructures. ADORBOND® CC PLUS Powder 10-30 µm can be used as a partial denture alloy as well as a metal-ceramic alloy.

Composition w_i

Co	%	63,6
Cr	%	24,8
W	%	5,5
Mo	%	5,0
Si	%	1,1

Properties

Particle size d _i	µm	10-30
Density ρ	g · cm ⁻³	8,2
Vickers hardness	HV 10	520
Linear thermal expansion coefficient α 25 - 500 °C	10 ⁻⁶ ·K ⁻¹	14,3
Linear thermal expansion coefficient α 20 - 600 °C	10 ⁻⁶ ·K ⁻¹	14,5
Melting range T _S - T _L	°C	1360-1420
Highest recommended firing temperature T _{Br,max}	°C	980
0,2-% Yield strength R _{p0,2}	MPa	1020
Modulus of elasticity E	GPa	275
Tensile elongation at break A ₅	%	5

Recommendations for Use

Digital Design Modelling

The design is carried out with suitable CAD software in consideration of dental rules. Ensure anatomically reduced framework shapes for subsequent ceramic veneering. Wall thickness after finishing: min. 0.3 mm. This ensures the clinically required stability. For implant prosthetic work with screw channels, the wall thickness should not be less than 0.5 mm. When setting the parameters in the CAD software, make sure that the finished object has a wall thickness of at least 0.44 mm. Wall thickness after finishing of the object: at least 0.3 mm. Design the connector as strong and high as possible (height: at least 3,5 mm, width: at least 2,5 mm; cross section min. 7 mm²). Ensure a sufficient connector cross section for pontics in the posterior region. Avoid undercuts and sharp edges. Particularly sensitive areas, such as at the transition from a connector to the terminal pontic, should be increased in wall thickness, at least 0.7 mm for orthodontic bands and even 1.5 mm for brackets.

SLM-Process

Avoid dust production while opening the powder container, transportation and filling into the SLM system! Use a suitable SLM system for the particle distribution of ADORBOND® CC PLUS Powder and carry out the SLM-Process according to the manufacturer's instructions for the SLM system. Non melted powder can be reused but should be sieved with a powder sieve (80 µm) first. Melted powder cannot be reused by e.g. casting.

Shot Peening

last the production platform with spherical ceramic material with Ø0,125-0,250 mm and 2,5-3,5 bar until a uniform saturation of the surface is visibly achieved.

Stress Relief Heat Treatment

After the laser sintering process has been carried out, the building panels must be heat-treated, which reduces internal stresses. A suitable furnace with inert gas (argon) or vacuum function is recommended for this purpose. The heat treatment needs to be carried out in a special furnace under argon atmosphere. Place the production platform in an inert gas box with an argon flow of 1-2 L/min and put the box into the cold furnace. Then:

- Within 60 min heat up to 450 °C and hold the temperature for 45 min.
- Within 45 min heat up to 750 °C and hold the temperature for 40-60 min
- Let the platform cool down. Open furnace door at max. 100 °C and remove the inert gas box and stop the argon flow.

Note: It is possible to carry out the test without inert gas, in which case increased attention must be paid to the fact that increased oxide formation may occur. The cooling phase should be carried out in the oven until room temperature is reached.

Firing of Ceramics

Use commercially available dental ceramics for cobalt based metal alloys with a suitable linear thermal expansion coefficient. Please follow the associated work instructions and cooling schemes given by the ceramic manufacturer. After the heat treatment and cooling down:

- Remove restorations from the plate using a band saw, rotating instruments or pliers. Use pliers to remove the rest of the supports. Carbide cutters are recommended for finishing of the object.
- Sand blast the surface by use of a pencil-blaster with aluminium oxide 100 µm or 250 µm at 3.5 bar.
- Ultrasonically clean the frame in distilled water or degrease with ethyl acetate.
- The oxide firing is optional, to be done at about 960 °C under vacuum for 5 minutes. Always remove the oxide layer after cleaning firing by sand blasting with aluminium oxide and degrease again.
Note: A clean surface is best to avoid bubbles in ceramics.
- The opaque is applied on the surface by a first thin wash firing and a second evenly covering opaque layer.
Before firing always let the opaque dry for 5-10 minutes at 600 °C.
- Firing and cooling should be carried out in accordance to the ceramic manufacturer's instructions.
- After every firing step (dentine bake, build-up and glazing) cooling phase until ca. 750 °C.

If you use a ceramic bonder, please observe the process steps of the respective manufacturer.

Finishing

After firing of the ceramic, polish the frame with suitable grinding and polishing instruments for dental alloys up to high gloss.

Soldering and Welding

Soldering before firing of the frame can be carried out with ADORBOND® CC-Lot and high temperature flux. The width of the solder gap should be 0,05-0,2 mm. After brazing, flux residues and metal oxides must be acidified and the surfaces cleaned by steam blasting. Attention! Work on implant interfaces must not be soldered! For welding with a laser use ADORBOND® CC-Laserdraht.

Safety Note

Metal dust is harmful to health and can cause irritation when inhaled and in contact with skin. When handling the powder, working out and blasting, use suction and wear protective goggles, protective gloves, respiratory mask (according to FFP3 - DIN EN 149) and protective clothing. After working with the powder or the restorations, thoroughly clean the skin contact areas, hands are particularly important.

Warranty

These application recommendations are based on own experiments and experiences and can therefore only be regarded as guidelines. The dentist or dental technician is responsible for the correct processing of this alloy.

Batch number

Refer to instructions for use

Manufacturer

Not for reuse